

Product Datasheet

Recombinant Human BAFFR Protein, His Tag (orb1550299)

Description	Recombinant Human BAFFR Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	BAFFR
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP2190, Recombinant Human Siglec-2/CD22 Protein is produced by mammalian expression system. The target protein is expressed with sequence (Ser7-Ala71) of human BAFFR (Accession #Q96RJ3-1) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.