

## Product Datasheet

# Recombinant Human BTLA Protein, Avi & His Tag (orb1550303)

<b>Description</b>	Recombinant Human BTLA Protein, Avi & His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	BTLA
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP2186, Recombinant Human PD-L1 Protein is produced by mammalian expression system. The target protein is expressed with sequence (Met1-Ser150) of human BTLA (Accession #Q7Z6A9) fused with a 6xHis Tag, Avi tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.