

Product Datasheet

Recombinant Human CD3 Delta Protein, His Tag (orb1550312)

Description	Recombinant Human CD3 Delta Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CD3 Delta
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP2177, Recombinant Human CD20 Protein is produced by mammalian expression system. The target protein is expressed with sequence (Phe22-Ala105) of human CD3 Delta (Accession #P04234) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.