

Product Datasheet

Recombinant Human CD30 Ligand/TNFSF8 Protein, Fc His Tag (orb1550319)

Description	Recombinant Human CD30 Ligand/TNFSF8 Protein, Fc His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CD30
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP2170, Recombinant Human IL-13 R alpha 1 Protein is produced by mammalian expression system. The target protein is expressed with sequence (Gln63-Asp234) of human CD30 Ligand/TNFSF8 (Accession #P32971-1) fused with an Fc at the N-terminus.
Expiration Date	6 months from date of receipt.