

Product Datasheet

Recombinant mouse FcRn Protein, Avi His Tag (orb1550322)

| | |
|------------------------|--|
| Description | Recombinant mouse FcRn Protein, Avi His Tag |
| Species/Host | Mammal |
| Reactivity | Mouse |
| Endotoxins | < 0.1 EU/ug by LAL method. |
| Target | FcRn |
| Note | For research use only |
| Purity | > 95% by SDS-PAGE. |
| Source | KMP2167, Recombinant Human Nectin-4 Protein is produced by mammalian expression system. The target protein is expressed with sequence (Ser22-Ser297 (FCGRT) & Ile21-Met119 (B2M)) of mouse FcRn (Accession #P01887&Q61559) fused with a 6xHis Tag, Avi tag at the C-terminus. |
| Expiration Date | 6 months from date of receipt. |