

Product Datasheet

Recombinant Human BAFF/TNFSF13B Protein, Fc His Tag (orb1550329)

Description	Recombinant Human BAFF/TNFSF13B Protein, Fc His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	BAFF
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP2160, Recombinant Human OX40/TNFRSF4 Protein is produced by mammalian expression system. The target protein is expressed with sequence (Ala141-Leu285 Trimer) of human BAFF/TNFSF13B (Accession #Q9Y275) fused with an Fc tag at the N-terminus.
Expiration Date	6 months from date of receipt.