

Product Datasheet

Recombinant Human TGF-beta RII / TGFBR2 Protein, mF & Avi Tag (orb1550339)

Description	Recombinant Human TGF-beta RII / TGFBR2 Protein, mF & Avi Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	TGFbeta RII
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2150, Recombinant Human TGF-beta RII / TGFBR2 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ile24-Asp159) of Human TGF-beta RII / TGFBR2 fused with a mFc His Tag and a Avi tag at C-terminal.
Expiration Date	6 months from date of receipt.