

Product Datasheet

Recombinant Human BCMA/TNFRSF17 Protein, Avi His Tag (orb1550340)

| | |
|------------------------|---|
| Description | Recombinant Human BCMA/TNFRSF17 Protein, Avi His Tag |
| Species/Host | Mammal |
| Reactivity | Human |
| Endotoxins | < 0.1 EU/ug by LAL method. |
| Target | BCMA |
| Note | For research use only |
| Purity | > 95% by SDS-PAGE. |
| Source | KMP2149, Recombinant Human 4-1BB Ligand/TNFSF9 Protein is produced by mammalian expression system. The target protein is expressed with sequence (Met1-Ala54) of human BCMA/TNFRSF17 (Accession #Q02223) fused with a 6xHis Tag, Avi tag at the C-terminus. |
| Expiration Date | 6 months from date of receipt. |