

Product Datasheet

Recombinant Human TNF alpha Protein, Avi His Tag (orb1550359)

Description	Recombinant Human TNF alpha Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	TNF alpha
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2130, Recombinant Human TNF alpha Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Val77-Leu233) of Human TNF alpha fused with His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.