

## Product Datasheet

# Recombinant Human TGF beta 2 Protein (orb1550362)

<b>Description</b>	Recombinant Human TGF beta 2 Protein
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	TGF beta 2
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
<b>Source</b>	KMP2127, Recombinant Human TGF beta 2 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Arg303-Ser414) of Human TGF beta 2 fused with no tag .
<b>Expiration Date</b>	6 months from date of receipt.