

Product Datasheet

Recombinant Human Transferrin R / CD71 Protein, Avi His Tag (orb1550368)

Description	Recombinant Human Transferrin R / CD71 Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Transferrin R
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2121, Recombinant Human Transferrin R / CD71 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Cys89-Phe760) of Human Transferrin R / CD71 fused with a His Tag and Avi tag at the N-terminal.
Expiration Date	6 months from date of receipt.