

Product Datasheet

Recombinant Human TIMP-2 Protein, His Tag (orb1550369)

Description	Recombinant Human TIMP-2 Protein, His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	TIMP2
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2120, Recombinant Human TIMP-2 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Cys27-Pro220) of Human TIMP-2 fused with a His Tag at the C- terminal.
Expiration Date	6 months from date of receipt.