

Product Datasheet

Recombinant Human VEGF165 Protein (orb1550376)

Description	Recombinant Human VEGF165 Protein
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	VEGF165
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2113, Recombinant Human VEGF165 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ala27-Arg191) of Human VEGF165 fused with no tag .
Expiration Date	6 months from date of receipt.