

Product Datasheet

Biotinylated Recombinant Human VEGF165 Protein, Biotinylated (orb1550377)

Description	Biotinylated Recombinant Human VEGF165 Protein, Biotinylated
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	VEGF165
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	KMP2112, Biotinylated Recombinant Human VEGF165 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ala27-Arg191) of Human VEGF165 fused with a His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.