

Product Datasheet

Recombinant Human P-Selectin/CD62P Protein, Avi His Tag (orb1550380)

Description	Recombinant Human P-Selectin/CD62P Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	PSelectin
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2109, Recombinant Human P-Selectin/CD62P Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Trp42-Ala771) of Human P-Selectin/CD62P fused with His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.