

Product Datasheet

Recombinant Human KDR / VEGFR2 Protein, mFc His Tag (orb1550383)

Description	Recombinant Human KDR / VEGFR2 Protein, mFc His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	KDR
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2106, Recombinant Human KDR / VEGFR2 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ala20-Glu764) of Human KDR / VEGFR2 fused with a mFc His Tag tag at the C-terminal.
Expiration Date	6 months from date of receipt.