

Product Datasheet

Recombinant Human CD70/CD27 Ligand Trimer Protein, His Tag (orb1550390)

Description	Recombinant Human CD70/CD27 Ligand Trimer Protein, His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CD70
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2099, Recombinant Human CD70/CD27 Ligand Trimer Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gln50-Pro193 Trimer) of Human CD70/CD27 Ligand Trimer fused with His Tag and Flag at the N-terminal.
Expiration Date	6 months from date of receipt.