

Product Datasheet

Recombinant Human siglec 9 Protein, Avi His Tag (orb1550392)

Description	Recombinant Human siglec 9 Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	siglec 9
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2097, Recombinant Human siglec 9 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gln18-Gly348) of Human siglec 9 fused with a His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.