

Product Datasheet

Biotinylated Recombinant Human siglec 10 Protein, hFAvi His Tag (orb1550394)

Description	Biotinylated Recombinant Human siglec 10 Protein, hFAvi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	siglec 10
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2095, Biotinylated Recombinant Human siglec 10 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Met17-Thr546) of Human siglec 10 fused with hFc and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.