

Product Datasheet

Biotinylated Recombinant Human ROR2 Protein, Avi His Tag (orb1550397)

Description	Biotinylated Recombinant Human ROR2 Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	ROR2
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2092, Biotinylated Recombinant Human ROR2 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Val34-Gly403) of Human ROR2 fused with His Tag and Avi tag at theC-terminal.
Expiration Date	6 months from date of receipt.