

## Product Datasheet

# Recombinant Human Siglec 10 Protein, hFAvi His Tag (orb1550400)

<b>Description</b>	Recombinant Human Siglec 10 Protein, hFAvi His Tag
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	Siglec 10
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by Tris-Bis PAGE
<b>Source</b>	KMP2089, Recombinant Human Siglec 10 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Met17-Thr546) of Human Siglec 10 fused with hFc and Avi tag at the C-terminal.
<b>Expiration Date</b>	6 months from date of receipt.