

Product Datasheet

Biotinylated Recombinant Human siglec 8 Protein, Avi His Tag (orb1550403)

Description	Biotinylated Recombinant Human siglec 8 Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	siglec 8
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2086, Biotinylated Recombinant Human siglec 8 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Met17-Ala363) of Human siglec 8 fused with a His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.