

## Product Datasheet

# Biotinylated Recombinant Human P-Selectin/CD62P Protein, Avi His Tag (orb1550407)

<b>Description</b>	Biotinylated Recombinant Human P-Selectin/CD62P Protein, Avi His Tag
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	PSelectin
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
<b>Source</b>	KMP2082, Biotinylated Recombinant Human P-Selectin/CD62P Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Trp42-Ala771) of Human P-Selectin/CD62P fused with His Tag and Avi tag at the C-terminal.
<b>Expiration Date</b>	6 months from date of receipt.