

Product Datasheet

Biotinylated Recombinant Human NKG2A Protein, Avi His Tag (orb1550408)

Description	Biotinylated Recombinant Human NKG2A Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	NKG2A
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	KMP2081, Biotinylated Recombinant Human NKG2A Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Arg100-Leu233) of Human NKG2A fused with His Tag and Avi tag at N-terminal.
Expiration Date	6 months from date of receipt.