

Product Datasheet

Recombinant Human siglec 6 Protein, Avi His Tag (orb1550413)

| | |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Recombinant Human siglec 6 Protein, Avi His Tag |
| Reactivity | Human |
| Endotoxins | < 1 EU/ug by LAL method. |
| Target | siglec 6 |
| Note | For research use only |
| Purity | > 95% by Tris-Bis PAGE, > 95% by SEC-HPLC |
| Source | KMP2076, Recombinant Human siglec 6 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gln27-Val331) of Human siglec 6 fused with a His Tag and Avi tag at the C-terminal. |
| Expiration Date | 6 months from date of receipt. |