

Product Datasheet

Biotinylated Recombinant Human PDGF R beta Protein, Avi His Tag (orb1550414)

Description	Biotinylated Recombinant Human PDGF R beta Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	PDGF R beta
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2075, Biotinylated Recombinant Human PDGF R beta Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Leu33-Phe530) of Human PDGF R beta fused with His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.