

Product Datasheet

Recombinant Human PD-L2/B7-DC Protein, Avi His Tag (orb1550415)

Description	Recombinant Human PD-L2/B7-DC Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	PDL2
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2074, Recombinant Human PD-L2/B7-DC Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Leu20-Thr220) of Human PD-L2/B7-DC fused with a His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.