

## Product Datasheet

# Recombinant Mouse PD-L1 Protein, hFc His Tag (orb1550417)

<b>Description</b>	Recombinant Mouse PD-L1 Protein, hFc His Tag
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	PDL1
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
<b>Source</b>	KMP2072, Recombinant Mouse PD-L1 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Phe19-Arg238) of Mouse PD-L1 fused with a hFc His Tag tag at the C-terminal.
<b>Expiration Date</b>	6 months from date of receipt.