

Product Datasheet

Recombinant Human NKG2A&CD94 Protein, Avi His Tag (orb1550419)

Description	Recombinant Human NKG2A&CD94 Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	NKG2A&CD94
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2070, Recombinant Human NKG2A&CD94 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (CD94(Ser34-Ile179)&NKG2A(Arg100-Leu233)) of Human NKG2A&CD94 fused with His Tag and Avi tag at the N-terminal.
Expiration Date	6 months from date of receipt.