

Product Datasheet

Recombinant Human NKp46 / NCR1 / CD335 Protein, mFc His Tag (orb1550433)

Description	Recombinant Human NKp46 / NCR1 / CD335 Protein, mFc His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	NKp46
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2056, Recombinant Human NKp46 / NCR1 / CD335 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gln22-Asn255) of Human NKp46 / NCR1 / CD335 fused with a mFc His Tag tag at the C-terminal.
Expiration Date	6 months from date of receipt.