

Product Datasheet

Biotinylated Recombinant Human Angiopoietin-2/ANGPT2 Protein, Avi His Tag (orb1550437)

Description	Biotinylated Recombinant Human Angiopoietin-2/ANGPT2 Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Angiopoietin2
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2052, Biotinylated Recombinant Human Angiopoietin-2/ANGPT2 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Lys275-Phe496) of Human Angiopoietin-2/ANGPT2 fused with a His Tag and Avi tag at the N-terminal.
Expiration Date	6 months from date of receipt.