

Product Datasheet

Recombinant Human Noggin Protein (orb1550440)

Description	Recombinant Human Noggin Protein
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Noggin
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	KMP2049, Recombinant Human Noggin Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gln28-Cys232) of Human Noggin fused with no tag .
Expiration Date	6 months from date of receipt.