

Product Datasheet

Recombinant Human Integrin alpha 5 beta 1 (ITGA5&ITGB1) Heterodimer Protein, Avi His Tag (orb1550443)

Description	Recombinant Human Integrin alpha 5 beta 1 (ITGA5&ITGB1) Heterodimer Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Integrin alpha 5 beta 1
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2046, Recombinant Human Integrin alpha 5 beta 1 (ITGA5&ITGB1) Heterodimer Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Phe42-Tyr995(P08648) acidic tail& Gln21-His Tag728(P05556) basic tail) of Human Integrin alpha 5 beta 1 (ITGA5&ITGB1) Heterodimer fused with His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.