

## Product Datasheet

# Biotinylated Recombinant Human MICB Protein, Avi His Tag (orb1550445)

<b>Description</b>	Biotinylated Recombinant Human MICB Protein, Avi His Tag
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	MICB
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
<b>Source</b>	KMP2044, Biotinylated Recombinant Human MICB Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ala23-Gly298) of Human MICB fused with His Tag and Avi tag at the C-terminal.
<b>Expiration Date</b>	6 months from date of receipt.