

Product Datasheet

Recombinant Human MICB Protein, Avi His Tag (orb1550447)

Description	Recombinant Human MICB Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	MICB
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2042, Recombinant Human MICB Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ala23-Gly298) of Human MICB fused with His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.