

Product Datasheet

Recombinant Human LILRB1 Protein, His Tag (orb1550452)

Description	Recombinant Human LILRB1 Protein, His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	LILRB1
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2037, Recombinant Human LILRB1 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Met1-Gly221) of Human LILRB1 fused with His Tag at the C- terminal.
Expiration Date	6 months from date of receipt.