

Product Datasheet

Recombinant Human IL-23 R Protein, mFc His Tag (orb1550459)

Description	Recombinant Human IL-23 R Protein, mFc His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IL23 R
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2030, Recombinant Human IL-23 R Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gly24-Gly355) of Human IL-23 R fused with mFc His Tag at the C-terminal.
Expiration Date	6 months from date of receipt.