

Product Datasheet

Recombinant Human IL-6R alpha Protein, Avi His Tag (orb1550468)

Description	Recombinant Human IL-6R alpha Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IL6R alpha
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 94% by SEC-HPLC
Source	KMP2021, Recombinant Human IL-6R alpha Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Leu20-Pro365) of Human IL-6R alpha fused with a His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.