

Product Datasheet

Recombinant Human IL-13 R alpha 2 Protein, hFc His Tag (orb1550474)

Description	Recombinant Human IL-13 R alpha 2 Protein, hFc His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IL13 R alpha 2
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2015, Recombinant Human IL-13 R alpha 2 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Asp27-Leu342) of Human IL-13 R alpha 2 fused with hFc His Tag tag at the C-terminal.
Expiration Date	6 months from date of receipt.