

## Product Datasheet

# Recombinant Human IL-17A&IL-17F heterodimer Protein, Avi His Tag (orb1550475)

<b>Description</b>	Recombinant Human IL-17A&IL-17F heterodimer Protein, Avi His Tag
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	IL17A
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
<b>Source</b>	KMP2014, Recombinant Human IL-17A&IL-17F heterodimer Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gly24-Ala155(IL-17A)&Arg31-Gln163) of Human IL-17A&IL-17F heterodimer fused with His Tag and Avi tag at the C-terminal.
<b>Expiration Date</b>	6 months from date of receipt.