

Product Datasheet

Recombinant Human IL7 Protein (orb1550476)

Description	Recombinant Human IL7 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IL7
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2013, Recombinant Human IL7 Protein is produced by E.coli expression system. The target protein is expressed with sequence (Asp26-His Tag177) of Human IL7 fused with no tag .
Expiration Date	6 months from date of receipt.