

Product Datasheet

Recombinant Human IL-27 Ra / TCCR Protein, hFc His Tag (orb1550477)

Description	Recombinant Human IL-27 Ra / TCCR Protein, hFc His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IL27 Ra
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	KMP2012, Recombinant Human IL-27 Ra / TCCR Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gly34-Lys516) of Human IL-27 Ra / TCCR fused with a hFc His Tag tag at the C-terminal.
Expiration Date	6 months from date of receipt.