

Product Datasheet

Biotinylated Recombinant Human IL-4 Protein, Avi His Tag (orb1550488)

Description	Biotinylated Recombinant Human IL-4 Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IL4
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	KMP2001, Biotinylated Recombinant Human IL-4 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Met1-Ser153) of Human IL-4 fused with a His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.