

Product Datasheet

Recombinant Human IL-2 Protein (orb1550494)

Description	Recombinant Human IL-2 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IL2
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1995, Recombinant Human IL-2 Protein is produced by E.coli expression system. The target protein is expressed with sequence (Ala21-Thr153) of Human IL-2 fused with no tag .
Expiration Date	6 months from date of receipt.