

Product Datasheet

Recombinant Human IL-17 RA / CD217 Protein, hFAvi His Tag (orb1550495)

Description	Recombinant Human IL-17 RA / CD217 Protein, hFAvi His Tag
Reactivity	Human
Endotoxins	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Target	IL17 RA
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	Recombinant Human IL-17 RA / CD217 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Leu33-Trp320) of Human IL-17 RA / CD217 fused with a hFc and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)