

## Product Datasheet

### Recombinant Human IL-18 Protein (orb1550500)

<b>Description</b>	Recombinant Human IL-18 Protein
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	IL18
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
<b>Source</b>	KMP1989, Recombinant Human IL-18 Protein is produced by E.coli expression system. The target protein is expressed with sequence (Tyr37-Asp193) of Human IL-18 fused with no Tag.
<b>Expiration Date</b>	6 months from date of receipt.