

## Product Datasheet

# Recombinant Human IL-3 Protein, Avi His Tag (orb1550504)

<b>Description</b>	Recombinant Human IL-3 Protein, Avi His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	IL3
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
<b>Source</b>	KMP1985, Recombinant Human IL-3 Protein is produced by E.coli expression system. The target protein is expressed with sequence (Ala20-Phe152) of Human IL-3 fused with a His Tag and Avi tag at the N-terminal.
<b>Expiration Date</b>	6 months from date of receipt.