

Product Datasheet

Biotinylated Recombinant Human IFN-gamma / IFNG Protein, Avi His Tag (orb1550506)

Description	Biotinylated Recombinant Human IFN-gamma / IFNG Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IFNgamma
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	KMP1983, Biotinylated Recombinant Human IFN-gamma / IFNG Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gln24-Gly161) of Human IFN-gamma / IFNG fused with a His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.