

Product Datasheet

Recombinant Human Guanylyl Cyclase C/GUCY2C Protein, Avi His Tag (orb1550517)

Description	Recombinant Human Guanylyl Cyclase C/GUCY2C Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Guanylyl Cyclase C
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1972, Recombinant Human Guanylyl Cyclase C/GUCY2C Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ser24-Gln430) of Human Guanylyl Cyclase C/GUCY2C fused with His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.