

Product Datasheet

Recombinant Human CD26/DPP4 Protein, hFc His Tag (orb1550522)

Description	Recombinant Human CD26/DPP4 Protein, hFc His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CD26
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1967, Recombinant Human CD26/DPP4 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Asn29-Pro766) of Human CD26/DPP4 fused with a hFc tag at the N-terminal.
Expiration Date	6 months from date of receipt.